



**Explore Genetic Technology to  
Improve the Quality of Life**



## **Pluslife FPV/ CPV Card**

Feline/Canine parvovirus

Rapid and accurate on-site POC molecular tests

# Current Challenges




Lab environment and skillful technician required

PCR testing



TAT=1-3 days

Independent clinical lab testing



Low sensitivity

Antigen testing

# Features




- FAST**  
TAT within **30** mins
- EASY**  
Hands-on time in **1~2** mins
- ACCURATE**  
Comparable to **PCR**
- ACCESSIBLE**  
Sample to result within **3** steps


## Pluslife FPV/CPV Card



## Test Procedure



Collect Sample



Mix Sample



Run Test

## Result Output



✓ Mini Dock



✓ PC Software



✓ Mobile App

# Application



Pet Hospital/Clinic





Pet Shop



Center for Animal Disease Control

# Specification




Pluslife Mini Dock	
Dimension	101 x 91 x 65 mm
Weight	210 g
Connection	Mobile phone or computer app
Operate temperature	15°C~28°C
Operate humidity	≤80%

Pluslife FPV/CPV Card	
Target pathogen	FPV/CPV
Sample type	Rectal swab
Turnaround time (TAT)	30 mins
Shelf life	13 months at 2°C~28°C

## Guangzhou Pluslife Biotech Co.,Ltd.

 Floor 3, Block E, Runhui Science & Technology Park, No.18 Shenzhou RD, Huangpu District, Guangzhou, Guangdong, China

 +86-20-31703986

 [market@pluslife.com](mailto:market@pluslife.com)

 [www.pluslife.com](http://www.pluslife.com)



Website



LinkedIn



Line